

Recombinant Mouse Galectin-1 (Lgals1)

Catalog No: #AP89677

Package Size: #AP89677-1 20ug #AP89677-2 100ug #AP89677-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Galectin-1 (Lgals1)
Brief Description	Recombinant Mouse Galectin-1 (Lgals1)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Target Name	His-tagged
Target Species	Mus musculus (Mouse)
Calculated MW	22.2 kDa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Target Sequence	Protein Length:Full Length of Mature Protein Expression Region: 2-135aa ACGLVASNLNLKPGECKVVRGEVASDAKSFVLNLGKDSNNLCLHFNPRFNAHGDANTIVCNTKEDGTWGTEH REPAFPFQPGSITEVCITFDQADLTIKLPDGHEFKFPNRLNMEAINYMAADGDFKIKCVAFE
Storage	" The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C."

Note: This product is for in vitro research use only and is not intended for use in humans or animals.