

# Human Ras-related protein Rap-2a (RAP2A) ELISA Kit

Catalog No: #EK20655

Package Size: #EK20655-1 48T #EK20655-2 96T



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Ras-related protein Rap-2a (RAP2A) ELISA Kit
Brief Description	Human Ras-related protein Rap-2a (RAP2A) ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.132ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.