

# DTX4 rabbit pAb

Catalog No: #49028



Package Size: #49028-1 50ul #49028-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DTX4 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of DTX4 at Human, Mouse
Immunogen Description	Synthesized peptide derived from human DTX4
Target Name	DTX4 KIAA0937 RNF155
Calculated MW	68kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.