

## Rabbit Anti-phospho-AMPK alpha-2 (Thr172) antibody

Catalog No: #C05078



Package Size: #C05078-1 50ul #C05078-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-phospho-AMPK alpha-2 (Thr172) antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB ELISA IHC IF
Species Reactivity	Hu Ms Rt
Target Name	phospho-AMPK alpha-2 (Thr172)
Other Names	phospho-AMPK alpha-2 (Thr172)
Accession No.	Uniprot:P54646Gene ID:5563
Calculated MW	64
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000

ELISA=1:5000-10000

IHC-P=1:100-500

IHC-F=1:100-500

IF=1:100-500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.