

Rat Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit

Catalog No: #EK20662

Package Size: #EK20662-1 48T #EK20662-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit
Brief Description	Rat Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.068ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.