

Recombinant Human Developmental pluripotency-associated protein 2 (DPPA2)

Catalog No: #AP89679

Package Size: #AP89679-1 20ug #AP89679-2 100ug #AP89679-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Developmental pluripotency-associated protein 2 (DPPA2)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 85% as determined by SDS-PAGE.
Target Sequence	<p>Expression Region:1-298aa</p> <p>Sequence Info:Full Length</p> <p>MSDANLDSSK KNFLEGEVDD EESVILTLVP VKDDANMEQM EPSVSSTSDV KLEKPKKYNP GHLLQTNEQF TAPQKARCKI PALPLPTILP PINKVCRDTL RDWCQQLGLS TNGKKIEVYL RLHRHAYPEQ RQDMPEMSQE TRLQRCSRKR KAVTKRARLQ RSYEMNERAE ETNTVEVITS APGAMLASWA RIAARAVQPK ALNCSIPVS VEAFMQASG VRWCVVHGRL LSADTKGWVR LQFHAGQAWV PTTHRRMISL FLLPACIFPS PGIEDNMLCP DCAKRNKMM KRLMTVEK</p>
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.