

Recombinant Human E3 ubiquitin/ISG15 ligase TRIM25 (TRIM25)



Catalog No: #AP89682

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89682-1 20ug #AP89682-2 100ug #AP89682-3 1mg

Description

Product Name	Recombinant Human E3 ubiquitin/ISG15 ligase TRIM25 (TRIM25)
Brief Description	Recombinant Human E3 ubiquitin/ISG15 ligase TRIM25 (TRIM25)
Purification	>85% (SDS-PAGE)
Target Name	The tag type will be determined during production process
Other Names	E3 ubiquitin/ISG15 ligase TRIM25; EFP; Estrogen responsive finger protein; Estrogen-responsive finger protein; RING finger protein 147; RNF 147; RNF147; TRI25; TRI25_HUMAN; TRIM 25; Trim25; Tripartite motif containing 25; Tripartite motif containing protein 25; Tripartite motif protein TRIM25; Tripartite motif-containing protein 25; Tripartite motifcontaining25; Ubiquitin/ISG15-conjugating enzyme TRIM25; Z147; Zinc finger protein 147 (estrogen responsive finger protein); Zinc finger protein 147; Zinc finger protein-147; ZNF 147; ZNF147
Target Species	Homo sapiens (Human)
Target Sequence	"Protein Length:Full length protein Expression Region: 1-630 MAELCPLAEE LSCSICLEPF KEPVTTPCGH NFCGSCNET WAVQGSPYLC PQCRAVYQAR PQLHKNTVLC NVVEQFLQAD LAREPPADVW TTPARASAPS PNAQVACDHC LKEAAVKTCL VCMASFQEH LQPHFDSPAF QDHPLQPPVR DLLRRKCSQH NRLREFFCPE HSECICHICL VEHKTCSPAS LSQASADLEA TLRHKLTVMY SQINGASRAL DDVRNRQQDV RMTANRKVEQ LQQEYTEMKA LLDASETST RKIKEEEKRV NSKFDTIYQI LLKKKSEIQT LKEEIEQSLT KRDEFEFLEK ASKLRGISTK PVYIPEVELN HKLIKGIHQ S TIDLKNEKQ CIGRLQEPTP SSGDPGEHDP ASTHKSTRPV KKSKEEKKS KKPPPVPALP SKLPTFGAPE QLVDLKQAGL EAAAKATSSH PNSTSLKAKV LETFLAKSRP ELLEYIYKVI LDYNTAHNKV ALSECYTVAS VAEMPQNYRP HPQRFTYCSQ VLGLHCYKKG IHYWEVELQK NNFCGVGICY GSMNRQGPES RLG SNSASWC VEFNTKISA WHNNVEKTL P STKATRVGV L NCNDHGFVIF FAVADKVHLM YKFRVDFTEA LYPAFWVFSA GATLSICSPK"
Storage	" The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C."

Note: This product is for in vitro research use only and is not intended for use in humans or animals.