

SLC35F1 Rabbit pAb

Catalog No: #49043



Package Size: #49043-1 50ul #49043-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC35F1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB, IHC-P, IHC-F, ICC, IF, ELISA
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human SLC35F1: 301-400/408
Other Names	C6orf169; dJ230I3.1; FLJ13018; S35F1_HUMAN; Slc35f1; Solute carrier family 35 member F1.
Calculated MW	45 kDa
Concentration	1mg/ml
Formulation	Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.