

Human Calsyntenin-1(CLSTN1) ELISA Kit

Catalog No: #EK20671



Package Size: #EK20671-1 48T #EK20671-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Calsyntenin-1(CLSTN1) ELISA Kit
Brief Description	Human Calsyntenin-1(CLSTN1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:35.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.