

Human Leucine-rich repeat transmembrane protein FLRT2 (FLRT2) ELISA Kit

Catalog No: #EK20683



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20683-1 48T #EK20683-2 96T

Description

Product Name	Human Leucine-rich repeat transmembrane protein FLRT2 (FLRT2) ELISA Kit
Brief Description	Human Leucine-rich repeat transmembrane protein FLRT2 (FLRT2) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL
Sensitivity:32.6pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.