

# Mouse Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit

Catalog No: #EK20669



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20669-1 48T #EK20669-2 96T

## Description

Product Name	Mouse Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit
Brief Description	Mouse Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:31.2-2000 pg/mL
Sensitivity:12.5pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.