

Human 17-ketosteroids (17-KS) ELISA Kit

Catalog No: #EK20696



Package Size: #EK20696-1 48T #EK20696-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human 17-ketosteroids (17-KS) ELISA Kit
Brief Description	Human 17-ketosteroids (17-KS) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:1-400 ng/mL

Sensitivity:0.5ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.