

Human Melanin concentrating hormone (MCH) ELISA Kit

Catalog No: #EK20701

Package Size: #EK20701-1 48T #EK20701-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Melanin concentrating hormone (MCH) ELISA Kit
Brief Description	Human Melanin concentrating hormone (MCH) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL
Sensitivity:0.113ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.