

## DNAJB9 Rabbit pAb

Catalog No: #49066



Package Size: #49066-1 50ul #49066-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNAJB9 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB, IHC-P, IHC-F, ICC, IF, ELISA
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human DNAJB9: 24-120/223
Other Names	DKFZP564F1862; DnaJ (Hsp40) homolog subfamily B member 9; DnaJ homolog subfamily B member 9; DNAJB9; DNJB9_HUMAN; Endoplasmic reticulum DnaJ homolog 4; ERdj4; MDG 1; Mdg-1; Microvascular endothelial differentiation gene 1; Microvascular endothelial differentiation gene 1 protein; MST049; MSTP049; OTTHUMP00000206840.
Calculated MW	23 kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.