

Glycogen synthase 2 antibody

Catalog No: #C05068



Package Size: #C05068-1 50ul #C05068-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycogen synthase 2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,IHC,IFoO ELISA
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Glycogen synthase 2
Target Name	Glycogen synthase 2
Other Names	Glycogen starch synthase liver; Glycogen starch synthase liver; EC 2.4.1.11; Glycogen synthase 2 liver; GYS2; GYS2_HUMAN.
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

WB=1:500-2000

IHC=1:100-500

IF=1:100-500

ELISA=1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.