

Rat Plasmin-antiplasmin complex (PAP) ELISA Kit

Catalog No: #EK20691



Package Size: #EK20691-1 48T #EK20691-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Plasmin-antiplasmin complex (PAP) ELISA Kit
Brief Description	Rat Plasmin-antiplasmin complex (PAP) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range: 3.12-200 ng/mL

Sensitivity: 1.39ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.