

PE anti-human CD140b

Catalog No: #248491



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | PE anti-human CD140b |
| Brief Description | PE anti-human CD140b |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 4A18A2 |
| Isotype | Mouse IgG1 |
| Applications | FC |
| Species Reactivity | Human |
| Other Names | Platelet-derived growth factor receptor, beta polypeptide, PDGFR1, PDGFRβ, PDGFRb, PDGF receptor bet |
| Formulation | 0.2%BSAε 0.1%Proclin 300η£-ι 'Eη • ηO ε ²ζ1AζA²oO PH7.2-7.4 |
| Storage | 4ε1 oO ι• Ωε δΩ• ε- |

Application Details

FC 5μl/ζ~'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.