

IRE1 Conjugated Rabbit pAb

Catalog No: #C55649

Package Size: #C55649-AF350 100ul #C55649-AF405 100ul #C55649-AF488 100ul #C55649-AF555 100ul #C55649-AF594 100ul #C55649-AF647 100ul #C55649-AF680 100ul #C55649-AF750 100ul #C55649-Biotin 100ul #C55649-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IRE1 Conjugated Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ERN1.
Other Names	ERN1;IRE1;IRE1P;IRE1a;hIRE1p
Accession No.	Swiss Prot:O75460GeneID:2081
Calculated MW	6kDa/109kDa
SDS-PAGE MW	130kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB1:500 - 1:2000

IHC1:50 - 1:200

IF/ICC1:50 - 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.