

# Canine Insulin (INS) ELISA Kit

Catalog No: #EK20712



Package Size: #EK20712-1 48T #EK20712-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Canine Insulin (INS) ELISA Kit
Brief Description	Canine Insulin (INS) ELISA Kit
Applications	ELISA
Species Reactivity	Canine

## Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.16ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.