

CHAC1 Polyclonal Antibody

Catalog No: #49244



Package Size: #49244-1 50ul #49244-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CHAC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	CHAC1 Polyclonal Antibody detects endogenous levels of protein.
Applications	WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	CHAC1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Target Name	CHAC1 BOTCH
Calculated MW	24kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

ELISA 1:5000-20000

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.