

Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit

Catalog No: #EK20729



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20729-1 48T #EK20729-2 96T

Description

| | |
|--------------------|---|
| Product Name | Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit |
| Brief Description | Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Mouse |

Application Details

Detect Range:12.5-800pg/mL

Sensitivity:5.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.