

Recombinant Human Polyadenylate-binding protein 2 (PABPN1)

Catalog No: #AP89684

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89684-1 20ug #AP89684-2 100ug #AP89684-3 1mg

Description

Product Name	Recombinant Human Polyadenylate-binding protein 2 (PABPN1)
Purification	>85% (SDS-PAGE)
Target Name	Tag type will be determined during the manufacturing process.
Other Names	Nuclear poly(A)-binding protein 1; OPMD; PAB2; PABII; PABP 2; pABP-2; PABP2; PABP2_HUMAN; PABPII; Pabpn1; poly(A) binding protein nuclear 1; Poly(A)-binding protein 2; Poly(A)-binding protein II; PolyA binding protein II; Polyadenylate-binding nuclear protein 1; Polyadenylate-binding protein 2
Target Species	Homo sapiens (Human)
Target Sequence	Protein Length: Full Length of Mature Protein Expression Region: 2-306 AAAAAAAAA AGAAGGRGSG PGRRRHLVPG AGGEAGEGAP GGAGDYGNGL ESEELEPEEL LLEPEPEPEP EEEPPRPRAP PGAPGPGPGS GAPGSQEEEE EPGLVEGDPG DGAIEDPELE AIKARVREME EEAELKELQ NEVEKQMNMS PPPGNAGPVI MSIEEKMEAD ARSIYVGNVD YGATAEELEA HFHGCGSVNR VTILCDKFSG HPKGFIYIEF SDKESVRTSL ALDESLFRGR QIKVIPKRTN RPISTTDRG FPRARYRART TNYNSSRSRF YSGFNSRPRG RVYRGRARAT SWYSPY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.