

PDZK1 rabbit pAb

Catalog No: #49639



Package Size: #49639-1 50ul #49639-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PDZK1 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of PDZK1 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human PDZK1
Target Name	PDZK1 CAP70 NHERF3 PDZD1
Calculated MW	57kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.