

Human Immunoglobulin Lambda Like Polypeptide 5 (IgLL5) ELISA Kit

Catalog No: #EK20737



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20737-1 48T #EK20737-2 96T

Description

Product Name	Human Immunoglobulin Lambda Like Polypeptide 5 (IgLL5) ELISA Kit
Brief Description	Human Immunoglobulin Lambda Like Polypeptide 5 (IgLL5) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL
Sensitivity:0.092ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.