

# Monkey Kidney Injury Molecule 1 (Kim1) ELISA Kit

Catalog No: #EK20738



Package Size: #EK20738-1 48T #EK20738-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Monkey Kidney Injury Molecule 1 (Kim1) ELISA Kit
Brief Description	Monkey Kidney Injury Molecule 1 (Kim1) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.