

Human Lysozyme g-like protein 1 (LYG1) ELISA Kit

Catalog No: #EK20740



Package Size: #EK20740-1 48T #EK20740-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Lysozyme g-like protein 1 (LYG1) ELISA Kit
Brief Description	Human Lysozyme g-like protein 1 (LYG1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.108ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.