

## Recombinant Human Erythropoietin protein(EPO)

Catalog No: #AP73384

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Erythropoietin protein(EPO)
Brief Description	Recombinant Protein
Host Species	Mammalian Cell
Purification	>98% as determined by SDS-PAGE and HPLC.
Immunogen Description	Expression Region:28-193aaSequence Info:Full Length of Mature Protein
Target Name	EPO
Accession No.	Uniprot ID: P01588
Target Species	Human
Calculated MW	21.0 kDa
Tag Info	Untagged
Target Sequence	APPRLICDSR VLERYLLEAK EAENITGCA EHCSLNENIT VPDTKVNIFYA WKRMEVGQQA VEVWQGLALL SEAVLRGQAL LVNSSQPWEP LQLHVDKAVS GLRSLTLLR ALGAQKEAIS PPDAASAAPL RTITADTFRK LFRVYSNFLR GKLKLYTGEA CRTGDR
Formulation	20mM Tris-HCl based buffer, pH7.6
Storage	Long Store -20°C, for extended storage, conserve at -20°C or -80°C. Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Application Details

Purity: Greater than 90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.