

# Beta-catenin(PT2250) mouse mAb

Catalog No: #61031

Package Size: #61031-1 50ul #61031-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-catenin(PT2250) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.