

CD43 mouse mAb(ABT027)

Catalog No: #61116

Package Size: #61116-1 50ul #61116-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD43 mouse mAb(ABT027)
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.