

Chromogranin A (ABT037) Mouse mAb

Catalog No: #61157

Package Size: #61157-1 50ul #61157-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chromogranin A (ABT037) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.