

CK13 mouse mAb(PT0076)

Catalog No: #61165

Package Size: #61165-1 50ul #61165-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CK13 mouse mAb(PT0076)
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.