

Cytokeratin 14 (PT1674) mouse mAb

Catalog No: #61223

Package Size: #61223-1 50ul #61223-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytokeratin 14 (PT1674) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.