

# Cytokeratin 19 (PT0087) mouse mAb

Catalog No: #61237

Package Size: #61237-1 50ul #61237-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                   |
|--------------------|-----------------------------------|
| Product Name       | Cytokeratin 19 (PT0087) mouse mAb |
| Host Species       | Mouse                             |
| Clonality          | Monoclonal                        |
| Isotype            | IgG                               |
| Applications       | IHC-p                             |
| Species Reactivity | Human                             |
| Storage            | Store at -20°C                    |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.