

# Cytokeratin 20 (CK20) mouse mAb(ABT413)

Catalog No: #61240

Package Size: #61240-1 50ul #61240-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytokeratin 20 (CK20) mouse mAb(ABT413)
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.