

Cytokeratin 4 (ABT-CK4) mouse mAb

Catalog No: #61243



Package Size: #61243-1 50ul #61243-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytokeratin 4 (ABT-CK4) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.