

Cytokeratin 6 (PT1886) mouse mAb

Catalog No: #61250

Package Size: #61250-1 50ul #61250-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------|
| Product Name | Cytokeratin 6 (PT1886) mouse mAb |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG |
| Applications | IHC-p |
| Species Reactivity | Human |
| Storage | Store at -20°C |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.