

# Fascin (PT1978) mouse mAb

Catalog No: #61288

Package Size: #61288-1 50ul #61288-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                           |
|--------------------|---------------------------|
| Product Name       | Fascin (PT1978) mouse mAb |
| Host Species       | Mouse                     |
| Clonality          | Monoclonal                |
| Isotype            | IgG                       |
| Applications       | IHC-p,IF/ICC              |
| Species Reactivity | Human                     |
| Storage            | Store at -20°C            |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.