

Gastrin mouse mAb(ABT174)

Catalog No: #61298

Package Size: #61298-1 50ul #61298-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------|
| Product Name | Gastrin mouse mAb(ABT174) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG |
| Applications | IHC-p |
| Species Reactivity | Human |
| Storage | Store at -20°C |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.