

HSP60 (PT1700) mouse mAb

Catalog No: #61338

Package Size: #61338-1 50ul #61338-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HSP60 (PT1700) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p,IF/ICC
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.