

Ki-67 mouse mAb(ABT001)

Catalog No: #61347

Package Size: #61347-1 50ul #61347-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ki-67 mouse mAb(ABT001)
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.