

p40 (DeltaNp63) mouse mAb(ABT020)

Catalog No: #61432

Package Size: #61432-1 50ul #61432-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	p40 (DeltaNp63) mouse mAb(ABT020)
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.