

# S100A1 mouse mAb(PT1428)

Catalog No: #61499

Package Size: #61499-1 50ul #61499-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                          |
|--------------------|--------------------------|
| Product Name       | S100A1 mouse mAb(PT1428) |
| Host Species       | Mouse                    |
| Clonality          | Monoclonal               |
| Isotype            | IgG                      |
| Applications       | IHC-p                    |
| Species Reactivity | Human                    |
| Storage            | Store at -20°C           |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.