

# TCL1 mouse mAb(ABT378)

Catalog No: #61526

Package Size: #61526-1 50ul #61526-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                        |
|--------------------|------------------------|
| Product Name       | TCL1 mouse mAb(ABT378) |
| Host Species       | Mouse                  |
| Clonality          | Monoclonal             |
| Isotype            | IgG                    |
| Applications       | IHC-p                  |
| Species Reactivity | Human                  |
| Storage            | Store at -20°C         |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.