

Tissue Inhibitor of Metalloproteinases 1(PT2236) mouse mAb

Catalog No: #61538

Package Size: #61538-1 50ul #61538-2 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tissue Inhibitor of Metalloproteinases 1(PT2236) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-p
Species Reactivity	Human
Storage	Store at -20°C

Note: This product is for in vitro research use only and is not intended for use in humans or animals.