

Human Putative pre-mRNA-splicing factor ATP-dependent RNA helicase DHX15 (DHX15) ELISA Kit



Catalog No: #EK20742

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20742-1 48T #EK20742-2 96T

Description

Product Name	Human Putative pre-mRNA-splicing factor ATP-dependent RNA helicase DHX15 (DHX15) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL

Sensitivity:0.125ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.