

# Human N-acetylglucosamine-1-phosphotransferase subunit gamma (GNPTG) ELISA Kit

Catalog No: #EK20759



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20759-1 48T #EK20759-2 96T

## Description

Product Name	Human N-acetylglucosamine-1-phosphotransferase subunit gamma (GNPTG) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.312-20 ng/mL
Sensitivity:0.103ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.