

# PE anti-mouse F4/80

Catalog No: #26848



Package Size: #26848-1 25T #26848-2 100T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PE anti-mouse F4/80
Isotype	rat IgG2b
Applications	FC
Species Reactivity	Mouse
Conjugates	PE
Other Names	EMR1, Ly71
Formulation	0.2%BSAε 0.1%Proclin 300η£·ι Έη • ηΟ ε ²ζΙΑζΑ²οΟ PH7.2-7.4
Storage	4 ε 'οΟ ι• Ωε δΩ• ε-

## Application Details

FC 5&mu;I/ζ~'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.