

Purified anti-mouse CD16/32

Catalog No: #AR27079



Package Size: #AR27079-1 100ug #AR27079-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse CD16/32
Isotype	rat IgG2b
Applications	FC
Species Reactivity	Ms
Other Names	FcγR III/II Receptor
Concentration	1.0mg/ml
Formulation	0.1 Proclin 300 phosphate buffered saline O PH7.2-7.4
Storage	Store at 2-8°C

Product Description

ε ζE©,12δE ζ

Note: This product is for in vitro research use only and is not intended for use in humans or animals.