



η» θ ζ θ ²ηΟ ε ²ζΑ²Cell Staining Buffer

Catalog No: #AR27080

Package Size: #AR27080 500ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	η» θ ζ θ ²ηΟ ε ²ζΑ²Cell Staining Buffer
Applications	FC
Formulation	PBS (pH7.4)+0.2%BSA+ 0.1% Proclin 300 οΟ 0.22 m ζ» ε ²θΩ ζ»
Storage	Store at 2-8°C

Application Details

η» θ ζ θ ²ηΟ ε ²ζΑ²οΟ Cell Staining

BufferοΟ η ²δΙ ζ• εΟ• η» θ ζξ ζ η» θ ζ θ ²θΩ η ²δΕ-η» θ η ζ ζΑ ε• η» θ ζ εΙθ° θ οΟ δΗ ε• —η ²δΙ ζ δ½ η η ²

δΕ»θΙ• ζ • ε PBSοΟ Proclin 300

οΟ BSAοΟ BSAδ½ δΕΙδΕ ηξ• θ η ½ε- ε ²δΙ ηΟ ε ²ζΑ²δΕ-ε• —δ» ε • ε° ζ δ½ η Ι• η ΗεΟ ζ ξη» ε• οΟ Proclin 300 ε• —δ» η»ζ • η» θ θ ²Ι• ζ ε η η³ε ζ ξΥ

Product Description

ε ζΕ©,12δΕ ζ

Note: This product is for in vitro research use only and is not intended for use in humans or animals.