

Monkey Amyloid Beta Peptide 1-40 (Ab1-40) ELISA Kit

Catalog No: #EK20747

Package Size: #EK20747-1 48T #EK20747-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monkey Amyloid Beta Peptide 1-40 (Ab1-40) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:125-8000pg/mL Sensitivity:43pg/mL Sample Type:serum, plasma, tissue homogenates, cell I

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.